



# Avoid Resampling Soil Vapors Confirm Tracer Gas in the Field Using a Leak Detector

By Irene DeGraff, Air Monitoring Product Marketing Manager, Russell Pellegrino, Director of Technical Services\*, and Kelli Steindl, GC Accessories Product Marketing Manager  
\* Centek Laboratories, LLC

- Confirm system integrity before sample collection.
- Minimize resampling by detecting leaks prior to sampling.
- Eliminate costly and time-consuming lab analysis of tracer gas.

Vapor intrusion occurs when pollutants from contaminated soil or ground water migrate into buildings and ambient air. Adverse health effects can result when vapors occur in high concentrations, or if toxic volatile organics are present. These compounds are monitored using a variety of sampling procedures, including soil vapor, sub-slab, indoor, and ambient air testing. Sample collection for volatile organic compounds (VOCs) typically is performed with an air canister and passive sampling kit according to EPA Method TO-15 or a similar method.

## Costly Detection in Lab Doesn't Prevent Resampling

The primary challenge in vapor intrusion monitoring is distinguishing vapor intrusion from other sources of exposure. In order to establish that VOCs are from soil vapor, rather than from the surrounding environment, sampling systems (ports) must be tested with tracer compounds, such as helium, and shown to be properly sealed. Sample collection system integrity can be demonstrated by including the tracer gas in the list of target analytes reported by the laboratory; however, if high levels are found the sample is rejected and costly resampling may result.

## Using a Leak Detector in the Field Saves Time and Money

Detection of tracer gas in the field is a cost-effective alternative to lab analysis that assures the integrity of the sampling system before sampling occurs. The Restek Leak Detector provides good screening of helium tracer gas at concentrations of 10%, the level at which sample port resealing is required. In addition, this unit is just a fraction of the cost of other field portable devices, such as photionization detectors, which may be too sensitive for screening purposes.


Real-time detection of helium tracer gas in the field using a Restek Leak Detector as shown in Figure 1 is a simple, inexpensive way to minimize resampling by establishing system integrity prior to sample collection. Centek Laboratories pioneered this technique and contributed to its inclusion in the New York State Department of Health method[1].

### References

1. New York State Department of Health, October 2006, Guidance for Evaluating Soil Vapor Intrusion in the State of New York, [http://www.nyhealth.gov/environmental/investigations/soil\\_gas/svi\\_guidance/docs/svi\\_main.pdf](http://www.nyhealth.gov/environmental/investigations/soil_gas/svi_guidance/docs/svi_main.pdf) (accessed August 27, 2010).

## Restek Mini-Cans

are ideal for both tracer gas transfer and introduction, as well as sample collection.



## Mini-Can Options

<b>Sizes</b>	400cc, 1000cc
<b>Valves</b>	Quick connect, diaphragm
<b>Interior Coating</b>	Electropolished, Siltek treated
<b>Sample Inlets</b>	Area, personal
<b>Flow ranges</b>	0.5-15 sccm

For a full product listing, visit [www.restek.com/air](http://www.restek.com/air)

**Prevent costly resampling—use a Restek Electronic Leak Detector to ensure sample collection system integrity prior to sampling.**

**Use a Restek Leak Detector to verify the sample port is sealed (Drawing A).**


1. Prepare sampling port by installing sample probe and shroud as described in NY DOH method.
2. Turn on Restek Electronic Leak Detector and allow it to equilibrate for a few seconds prior to use.
3. Insert leak detector probe tip into the open end of the tubing connected to the sealed sample port.
4. Wait 10 seconds and inject a charge of helium into the open space of the shroud.
5. Wait several minutes. An alarm will sound if helium is detected at >10%, indicating a leak in the sample port.





Photos courtesy of Centek Laboratories, LLC

**Collect sample using any Restek air canister or mini-can (Drawing B).**

**Drawing A: Tracer Gas Detection**



**Drawing B: Sample Collection**



**NEW!**

# Restek Electronic Leak Detector

**Protect your instrument and analytical column!**

High temperature GC methods are extremely sensitive to carrier gas impurities such as water and oxygen. Make sure you have clean carrier gas and frequently check connections and injection system fittings for leaks using a Restek Electronic Leak Detector.

Leak Detector with Universal Charger Set  
(US, European, Australian)

Cat.# 22839 . . . . . \$899